Automatic Seed Counter

Brand **Master Science** | Unit **Pc**

Automatic Seed Counter

MSASCCPLUSSWISS Automatic Seed Counter:Specification: Model: MSASCCPLUSSWISS Count Range: 1 to 12 mm Count Accuracy: ±2‰ Count Speed: ≥1000 seeds/3 minutes (for Sesame) Counting Range: 1 to 99,999 seeds Noise Level: <70 dB Power Supply: AC 220/110V, 50/60Hz Operation Environment: -10°C to 50°C, <85% RH External Size (W × D × H): 285 × 210 × 176 mm Packing Size (W × D × H): 430 × 280 × 230 mm Gross Weight: 9.5 kg

Request a Quotation
Contact us for pricing, configuration options and delivery details.

support@themsci.com
www.themsci.com

Scan to view online

Technical Specifications

PARAMETER	SPECIFICATION
Count Range	1 to 12 mm
Count Accuracy	±2‰
Count Speed	≥1000 seeds/3 minutes (for Sesame)
Counting Range	1 to 99,999 seeds
Noise Level	<70 dB
Power Supply	AC 220/110V, 50/60Hz
Operation Environment	-10°C to 50°C, <85% RH
External Size (W × D × H)	285 × 210 × 176 mm
Packing Size (W × D × H)	430 × 280 × 230 mm
Gross Weight	9.5 kg