Autosampler Vials 2 ml - HPLC GC Sample Vials - 12 × 32 mm - 9-425 Screw Thread

Model **MSC0000008SWISS, MSC0000013SWISS** | Brand **Master Science**








Master Science 2 ml Autosampler Vials are glass sample vials for liquid and gas chromatography. The 12 × 32 mm body with a wide-opening 9-425 screw thread is the format accepted by common autosamplers.

KEY HIGHLIGHTS

- Easy filling
- Protection from light
- Stable in the rack
- Quick volume check

Request a Quotation
Contact us for pricing, configuration options and delivery details.

support@themsci.com
www.themsci.com

Scan to view online

Overview

They are supplied in clear glass and in amber glass for light-sensitive samples, both with a flat bottom and graduation marks. Each pack contains 100 vials.

Key Benefits

Fits standard trays
Master Science 2 ml Autosampler Vials use the 12 × 32 mm format.

Easy filling
wide 9-425 opening.

Protection from light
amber option.

Stable in the rack
flat bottom.

Quick volume check
printed graduations.

They suit the autosampler of the Master Science rapid separation HPLC system.

Technical Specifications

GENERAL	
Brand	Master Science
Model Number	MSC0000008SWISS, MSC0000013SWISS
SKU	MSC0000008SWISS, MSC0000013SWISS
Brand	Master Science
Model Number	MSC0000008SWISS (clear), MSC0000013SWISS (amber)
SKU	MSC0000008SWISS, MSC0000013SWISS
Volume	2 ml
Material	Glass
Neck Type	Screw thread, 9-425 wide opening
Colour	Clear / amber
Bottom	Flat
Application	LC and GC
Package	100 pcs per pack

DIMENSIONS & WEIGHT	
Vial Size	12 × 32 mm

HOW MASTER SCIENCE 2 ML AUTOSAMPLER VIALS SUPPORT ROUTINE INJECTIONS	
Vials that do not fit the tray	standard 12 × 32 mm size.
Light-sensitive analytes degrading	amber glass.
Narrow necks that are hard to fill	wide opening.
Running out mid-batch	100 vials per pack.

Applications

Master Science 2 ml Autosampler Vials are used for
HPLC sample injection GC sample injection General sample storage in chromatography laboratories

Design

- Wide-opening 9-425 screw neck
- 12 × 32 mm outer dimensions
- Flat base

Choice

- Clear glass
- Amber glass