Dry Bath

Brand **Master Science** | Unit **Pc**

Dry Bath

Model: MSDBI1001SWISS Standard: 2 blocks Temperature Range: RT.+5°C ~ 160°C Temperature Accuracy: ±1°C (@120°C) Temperature Uniformity: ±0.5°C Display Accuracy: 0.1°C Time Range: 1-99h59min / ∞ Heating Time: ≤15 min (from 25°C to 160°C) Max Power: 500W Voltage: AC 220V, 50Hz Size: 260 × 220 × 100 mm Package Size: 360 × 310 × 230 mm Gross Weight: 6 kg

Request a Quotation
Contact us for pricing, configuration options and delivery details.

support@themsci.com
www.themsci.com

Scan to view online

Technical Specifications

PARAMETER	SPECIFICATION
Standard	2 blocks
Temperature Range	RT.+5°C ~ 160°C
Temperature Accuracy	±1°C (@120°C)
Temperature Uniformity	±0.5°C
Display Accuracy	0.1°C
Time Range	1-99h59min / ∞
Heating Time	≤15 min (from 25°C to 160°C)
Max Power	500W
Voltage	AC 220V, 50Hz
Size	260 × 220 × 100 mm

Standard Accessories Included

- 360 × 310 × 230 mm
- Gross Weight: 6 kg