
Electric Balance, 20Kg

Model **MSBE20001SWISS** | Brand **Master Science** | Unit **Pc**



Electric Balance, 20Kg (Model: MSBE20001SWISS): Technical Specifications: • Capacity: 0 ~ 20,000 g (20 kg) • Readability: 1 g • Repeatability: 1 g • Linear Error: 1 g • Scale Size: 255 x 190 mm • Gross Weight: 4.5 kg • Package Size: 460 x 395 x 250 mm Key Features: 1. High Capacity: Measures up to 20 kg, making it ideal for various industrial, laboratory, and commercial applications. 2.

KEY HIGHLIGHTS

- High Capacity

Request a Quotation
Contact us for pricing, configuration options and delivery details.

support@themsci.com
www.themsci.com

Scan to view online

Key Features

High Capacity
Measures up to 20 kg, making it ideal for various industrial, laboratory, and commercial applications.

The scale size of 255 x 190 mm allows for easy placement and use.

This electric balance is designed for precise weight measurements in laboratories, production facilities, and research environments.

Technical Specifications

PARAMETER	SPECIFICATION
Capacity	0 ~ 20,000 g (20 kg)
Readability	1 g
Repeatability	1 g
Linear Error	1 g
Scale Size	255 x 190 mm
Gross Weight	4.5 kg
Precision	Offers a high level of accuracy with 1 g readability, repeatability, and linear error.
Durability	Robust build with a manageable weight of 4.5 kg, suitable for extended use.

Standard Accessories Included

- 460 x 395 x 250 mm