# Glass Bead Sterilizer - Digital Dry Instrument Sterilizer - 100-300°C - 150g Bead Capacity

Model **MSBGS-300LSWISS** | Brand **Master Science**






The Master Science Glass Bead Sterilizer is a small benchtop unit for quick sterilization of forceps, scissors, needles and other fine surgical or dissecting instruments. Its chamber is filled with glass beads that are heated to as much as 300°C, and the tips of the instruments are pushed into the hot beads for about 10 to 15 seconds.

KEY HIGHLIGHTS

- Seconds per instrument
- No chemicals or flames
- Hood-friendly
- Saves bench space

AT A GLANCE

| | |
|---|---|
| **100 to 300°C**<br>TEMPERATURE RANGE | **1 min to 99 h 59 min, or continuous**<br>TIMER RANGE |
| **40 × 80 mm**<br>CHAMBER SIZE (Ø × DEPTH) | **±5°C**<br>TEMPERATURE ACCURACY |

**Request a Quotation**
Contact us for pricing, configuration options and delivery details.

support@themsci.com
www.themsci.com

Scan to view online

## Overview

Because it works with dry heat only, there are no gases, fumes, open flames or liquid disinfectants, so the Master Science Glass Bead Sterilizer can be used safely inside a laminar flow hood. A digital controller with LCD shows the real-time temperature and countdown, and over-temperature protection is built in.

Between-sample sterilization of small tools is often slow or unsafe. The Master Science Glass Bead Sterilizer helps with:

## Key Benefits

**Seconds per instrument**
Small instruments are sterilized in roughly 10 to 15 seconds once the beads are hot.

**No chemicals or flames**
Dry heat avoids gases, fumes and liquid disinfectants.

**Hood-friendly**
Safe for use in laminar flow hoods where open flames are not wanted.

**Saves bench space**
A compact footprint of 176 × 135 mm fits beside the work area.

For sterilizing inoculating loops, see our bacti-cinerator sterilizer.

## Technical Specifications

| GENERAL | |
|---|---|
| Brand | Master Science |
| Model Number | MSBGS-300LSWISS |
| SKU | MSBGS-300LSWISS |
| Brand | Master Science |
| Model | MSBGS-300LSWISS |
| SKU | MSBGS-300LSWISS |
| Temperature Range | **100 to 300°C** |
| Temperature Accuracy | **±5°C** |
| Timer Range | **1 min to 99 h 59 min, or continuous** |
| Heating Time | 25 min (room temperature to 300°C) |
| Glass Bead Capacity | **150 g** |
| Lid Material | Stainless steel |
| General | Max. Power: 120 W |
| Display Resolution | 1°C |
| Voltage | AC 110 V or 220 V, 50/60 Hz |

| DIMENSIONS & WEIGHT | |
|---|---|
| Chamber Size (Ø × Depth) | 40 × 80 mm |
| Dimensions | 176 × 135 × 189.5 mm |

| PROBLEMS THE MASTER SCIENCE GLASS BEAD STERILIZER SOLVES | |
|---|---|
| Slow turnaround between samples | Hot beads sterilize instrument tips in seconds. |
| Open flames in the hood | Dry bead heating needs no flame or alcohol. |
| Chemical residues | No liquid disinfectants touch the instruments. |
| Uncertain temperature | The digital display shows the actual bead temperature at all times. |

## Applications

**The Master Science Glass Bead Sterilizer is used for**

Sterilizing forceps, scissors and needles between samples Dissection and tissue culture work Plant and animal cell culture labs Dental and veterinary small instruments Microbiology work in laminar flow hoods

## Control and Display

- Digital temperature control up to 300°C
- LCD with real-time temperature display
- Countdown timer display
- Over-temperature protection

## Construction

- Stainless steel inner container
- Stainless steel lid
- Compact benchtop housing

## Standard Accessories Included

- Glass bead sterilizer main unit
- 1 bag of glass beads
- 3 fuses