
Refrigerated Centrifuge, Angle Rotor 12×15 ml, Speed 6000 rpm

Brand **Master Science** | Unit **Pc**



This centrifuge is equipped with a high-capacity rotor and a refrigeration system, ideal for applications requiring precise temperature control and speed accuracy.

Request a Quotation
Contact us for pricing, configuration options and delivery details.

support@themsci.com
www.themsci.com

Scan to view online

Technical Specifications

PARAMETER	SPECIFICATION
RotorAngle Rotor 12×10/15 ml	
Max. Speed6000 rpm	
Speed Accuracy±10 rpm	
Max. RCF (Relative Centrifugal Force)5290×g	
Max. Volume4×500 ml	
Refrigeration SystemCFC-free compressor, R404a refrigerant	
Timing Range1 second to 99 hours 59 seconds / Continuous / Short-span	
Temperature Range-20℃ to +40℃	
Temperature Accuracy±1℃	
Noise Level≤65 dB	
Power SupplyAC 220V, 50/60Hz	
Rated Power1200W	
External Dimensions768×672×444 mm	