Ultrasonic Bath /Ultrasonic Cleaner, 30L, 80 degree C

Model **MSVTD30L80DCNSWISS** | Brand **Master Science** | Unit **Pc**





Ultrasonic Cleaner MSVTD30L80DCNSWISS– Extreme Technology Features & Functions: Capacity: 30 liters, suitable for cleaning large items efficiently. Ultrasonic Frequency: Operates at 40KHz for effective ultrasonic cleaning. Adjustable Power & Time: Power range from 10% to 100%, time setting adjustable from 1 to 99 minutes. Large LCD Display: Displays precise time, temperature, and power levels for easy control.

KEY HIGHLIGHTS

- Large LCD Display
- Sweep Function
- Memory Function
- Energy-Saving Mode
- Low Noise Operation

Request a Quotation
Contact us for pricing, configuration options and delivery details.

support@themsci.com
www.themsci.com

Scan to view online

Key Features

Large LCD Display
Displays precise time, temperature, and power levels for easy control.

Sweep Function
Ensures uniform cleaning action by distributing ultrasonic energy evenly.

Memory Function
Remembers the last settings for convenient future use.

Energy-Saving Mode
Sleeping mode helps conserve energy.

Low Noise Operation
Thanks to built-in sound insulation, the machine runs quietly.

Includes a carrying handle and drain valve for easy usage and cleaning.

Material & Build
Top-Quality Stainless Steel: Durable stainless steel basket, body, tank, and cover for long-lasting use.

Flow Control Valve
Convenient drain valve for easy rinsing after cleaning.

Cleaning Basket
Comes with a basket for handling delicate items.

Technical Specifications

PARAMETER	SPECIFICATION
Capacity30 liters, suitable for cleaning large items efficiently.	
Ultrasonic FrequencyOperates at 40KHz for effective ultrasonic cleaning.	
Ultrasonic PowerMax 480W, adjustable in 10% steps.	
Heating PowerMax 500W.	
Capacity30 liters.	
Tank Dimensions500mm x 300mm x 150mm.	
External Dimensions530mm x 330mm x 340mm.	
Number of Transducers8.	
Drain ValveYes.	
Carrying HandleYes.	
Ultrasonic Frequency40KHz.	
Temperature SettingUp to 80°C	
Timer SettingAdjustable from 1 to 99 minutes.	
Adjustable Power & TimePower range from 10% to 100%, time setting adjustable from 1 to 99 minutes.	
Heating & Degas FunctionsBuilt-in heating up to 80°C, degas function for better cleaning performance.	

Standard Accessories Included

- Cleaning basket.
- This Ultrasonic Cleaner is ideal for professional cleaning tasks requiring precision, power, and versatility, combining efficient ultrasonic technology with user-friendly features for high-quality results.