Ultrasonic Bath /Ultrasonic Cleaner, 22L

Model **MSVTD3351NiSWISS** | Brand **Master Science** | Unit **Pc**





Ultrasonic Cleaner MSVTD3351NSWISS – Extreme Technology Features & Functions: Capacity: 22 liters, suitable for cleaning large items efficiently. Ultrasonic Frequency: Operates at 40KHz for effective ultrasonic cleaning. Adjustable Power & Time: Power range from 10% to 100%, time setting adjustable from 1 to 99 minutes. Large LCD Display: Displays precise time, temperature, and power levels for easy control.

KEY HIGHLIGHTS

- Large LCD Display
- Heating & Degas Functions
- Sweep Function
- Memory Function
- Energy-Saving Mode

Request a Quotation
Contact us for pricing, configuration options and delivery details.

support@themsci.com
www.themsci.com

Scan to view online

Key Features

Large LCD Display
Displays precise time, temperature, and power levels for easy control.

Heating & Degas Functions
Built-in heating up to 60°C (option for 80°C), degas function for better cleaning performance.

Sweep Function
Ensures uniform cleaning action by distributing ultrasonic energy evenly.

Memory Function
Remembers the last settings for convenient future use.

Energy-Saving Mode
Sleeping mode helps conserve energy.

Low Noise Operation
Thanks to built-in sound insulation, the machine runs quietly.

Includes a carrying handle and drain valve for easy usage and cleaning.

Flow Control Valve
Convenient drain valve for easy rinsing after cleaning.

Cleaning Basket
Comes with a basket for handling delicate items.

Technical Specifications

GENERAL	
General	Capacity22 liters, suitable for cleaning large items efficiently.
General	Ultrasonic FrequencyOperates at 40KHz for effective ultrasonic cleaning.
General	Ultrasonic PowerMax 480W, adjustable in 10% steps.
General	Heating PowerMax 500W.
General	Capacity22 liters.
General	Tank Dimensions500mm x 300mm x 150mm.
General	External Dimensions530mm x 330mm x 340mm.
General	Number of Transducers8.
General	Drain ValveYes.
General	Carrying HandleYes.
General	Ultrasonic Frequency40KHz.
General	Temperature SettingDefault 60°C
General	Timer SettingAdjustable from 1 to 99 minutes.
General	Adjustable Power & TimePower range from 10% to 100%, time setting adjustable from 1 to 99 minutes.

MATERIAL & BUILD	
Top-Quality Stainless Steel	Durable stainless steel basket, body, tank, and cover for long-lasting use.

Standard Accessories Included

- Cleaning basket.
- This Ultrasonic Cleaner is ideal for professional cleaning tasks requiring precision, power, and versatility, combining efficient ultrasonic technology with user-friendly features for high-quality results.